Incudigit-TFT

2001261	19L
2001262	36L
2001263	52L
2001264	80L
2001265	150L

Digitheat-TFT

2001251	19L
2001252	36L
2001253	52L
2001254	80L
2001255	150L

Digitronic-TFT

2005163	33L
2005165	47L
2005167	76L
2005169	145L

(con puerta vidrio)

2005164	33L
2005166	47L
2005168	76L
2005170	145L

Incubig-TFT

2000238	288L	2000239	400L	2000240	720L
---------	------	---------	------	---------	------

ESTUFAS DE CULTIVO / ESTUFAS DESECACIÓN
CULTURE OVENS / DESICCATION OVENS

Indice

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Index

1. General Information.....	14
2. Packing list	14
3. Technical specification	14
4. Information about safety in the use of the equipment.....	16
4.1 Risks the operator is subjected to	16
4.2 Staff qualification.....	16
4.3 Use with hazardous substances	16
4.4 Maintenance.....	16
4.5 Warranty.....	16
5. Equipment Description	17
6. Installation	17
6.1 Location.....	17
6.2 Connection to the mains	17
7. Operation.....	17
7.1 Fast temperature setting and start	18
7.2 Setting screen	19
7.3 TASK screen	21
7.4 PROGRAMMING screen	22

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1. General Information

- Handle the parcel with care. Unpack and check that the contents coincide with the packing list. If any part is damaged or missing, please advise the distributors immediately.
- Do not install or use the equipment without reading this handbook first. This instruction manual must be attached to the equipment and must be available for all users.
- If you have any doubts or enquiries, contact with J.P. SELECTA's technical service.
- Do not use the equipment in explosive, flammable or corrosive atmospheres.

2. Packing list

Model:	2001261	2001262	2001263	2001264	2001265
Shelves guides:	2x 2000011	2x 2000012	2x 2000012	2x 2000013	2x 2000015
Shelves:	2x 2000021	2x 2000022	2x 2000024	2x 2000023	2x 2000025

Model:	2000238	2000239	2000240
Shelves guides:	2x 2002371	--	--
Shelves:	2x 2002372	2x 2000063	2x 2000064

Model:	2001251	200125	2001253	2001254	2001255
Shelves guides:	2x 2000011	2x 2000012	2x 2000012	2x 2000013	2x 2000015
Shelves:	2x 2000021	2x 2000022	2x 2000024	2x 2000023	2x 2000025

Model:	2005163	2005165	2005167	2001269
	(glass door) 2005164	2005166	2005168	2005170
Shelves guides:	2x 2000012	2x 2000033	2x 2000013	2x 2000015
Shelves:	2x 2000072	2x 2000073	2x 2000074	2x 2000075

3. Technical specification

Type:	Culture oven				
Model:	2001261	2001262	2001263	2001264	2001265
Chamber space: (L)	19	36	52	80	150
Interior: High/wide/depth (cm)	30/25/25	40/30/30	33/47/33	50/40/40	50/60/50
External : High/wide/depth (cm)	51/57/49	60/62/54	53/79/57	70/72/64	70/92/74
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	170	225	250	300	525
Weigth: (Kg)	26	36	46	54	75
Fuse:	Replaceable (Type: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 80°C				
Display resolution:	0.1°C				
Timed operation:	week days, hh.mm				
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				

Type:	Culture oven		
Model:	2000238	2000239	2000240
Chamber space: (L)	288	400	720
Interior: High/wide/depth (mm)	800/600/600	1000/800/500	1200/1000/600
External: High/wide/depth (mm)	970/910/760	1300/1140/750	1520/1840/850
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)		
Electrical power (W):	570	1100	1600
Weight: (Kg)	87	160	225
Fuse:	Reemplazable (Type: 10x5mm)		
Temperature range:	From ambient + 5°C ... up to ... 80°C		
Display resolution:	0.1°C		
Timed operation:	Week, hour, minute		
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.		
Data recording:	USB memory device		

Type:	Drying stove (natural convection)				
Model:	2001251	2001252	2001253	2001254	2001255
Chamber space: (L)	19	36	52	80	150
Interior: High/wide/depth (cm)	30/25/25	40/30/30	33/47/33	50/40/40	50/60/50
External: High/wide/depth (cm)	50/60/44	60/65/49	53/82/52	70/74/59	70/95/68
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	600	900	1000	1200	2100
Weight: (Kg)	24	35	44	59	73
Fuse:	Replaceable (Type: 10x5mm)				
Temperature range:	From ambient + 5°C ... up to ... 250°C				
Display resolution:	1°C				
Timed operation:	Week days, hh.mm				
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				

Type:	Drying stove (forced convection)				
Modelo:	2005163	2005165	2005167	2001269	
	(glass door) 2005164	2005166	2005168	2005170	
Chamber space: (L)	33	47	76	145	
Interior: High/wide/depth (cm)	40/28/30	33/45/32	50/38/40	50/58/50	
External: High/wide/depth (cm)	60/65/55	53/81/58	70/75/65	70/95/72	
Power supply:	230 VAC 50Hz/60Hz (according to characteristics label)				
Electrical power (W):	1200	1200	1600	2000	
Weight (Kg):	40	50	64	79	
Fuse:	Reemplazable (Tipo: 10x5mm)				
Temperature range:	Desde... Tª ambiente + 5°C ... hasta ... 250°C				
Display resolution:	1°C				
Timed operation:	Week, hour, minute				
Temperature Profiles:	Up to 10 profiles 6 segments. Segment from 1minute to 99hour.				
Data recording:	USB memory device				

4. Information about safety in the use of the equipment

Hazard risk

Electric risk

Risk of high temperature contact

4.1 Risks the operator is subjected to

- Risk of burns
- Electrical risk

4.2 Staff qualification

This equipment may only be used by properly trained personnel who have read and understood these instructions or have been properly trained in this equipment operation.

4.3 Use with hazardous substances

Do not treat flammable or explosive materials or those which chemically react with violence.

ATTENTION!

Burning risk

Take proper precautions
(use of thermal gloves, clamps,...)
to extract samples from inside the oven.

It can be at a high temperature.

4.4 Maintenance

Unplug the equipment from the mains before removing the oven's cover to handle inside.

The manipulation of the equipment internal electronic circuits by unauthorized personnel can cause irreparable damage. Take it to one of J.P.SELECTA's authorized technical services.

It is recommended to clean the accessible parts of the oven with a cloth and a cleaner, or a similar solution. Do not use aggressive solvents.

4.5 Warranty

This product is guaranteed for one year. It does not cover damages caused by an incorrect use or causes beyond the control of J. P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

Fig 1. Temperature profile

5. Equipment Description

Ovens are used in sample thermal treatments, for its desiccation, cured and culture incubation.

Giving the heat treatment possibility by following a "temperature profile" allows avoiding cracks, breaks, color changes, etc.

A temperature profile is a succession of "ramps" and "dells."

A ramp is specified by a slope of °C/minute or °C/hour.

A dell is specified by a temperature and a time.

The device has a temperature log in an USB memory drive.

6. Installation

Unpack the equipment and visually check that the contents are in good conditions. Keep the package for some days.

6.1 Location

Place the equipment on a flat level stable surface, suitable for the equipment weight.

To work in a comfortable way, please leave a free space of 50cm around the equipment and a minimum of 20cm at the rear part.

6.2 Connection to the mains

Plug the equipment to a near socket suitable for the equipment power. Use the connection cable supplied or a similar one.

For your safety, the equipment must be earthed.

Before plugging the machine, check that the voltage supply and the frequency corresponds to the one indicated in the equipment characteristics plate.

7. Operation

All the equipment functions are set with the touch screen. Gently press with your finger in the center of each icon.

When starting the computer, the starting screen will be shown (Fig 7.1).

Fig 7.1. Initial screen

- **SETTINGS:**

- Time and date: menu to set the date and time.
- PID and Autotuning: menu to set the parameters of temperature control.
- Downloading recorded data.
- Temperature correction.
- Temperature alarm limits.

- **PROGRAMMING**

- Allows you to program the temperature profile, by defining temperatures, times and ramps.


- **TASK**

- Allows you to select and launch the operating temperature or a temperature profile.


7.1 Quick start to work at a fixed temperature


1 Press on «OPERATION».

2 Press over temperature panel. The setpoint indicator will change to blue.

3 Press   to select the setpoint temperature. Wait a few seconds to confirm the new value.

4 Press on  to start.


5 Press on  to end the cycle.

NOTE: at the end of the cycle, a message appears saying that the cycle has been completed by the user. To delete it, click on the message.

NOTE: If the oven is turned off with the ON/OFF switch and without pressing STOP at the end of the cycle, a power failure error message will appear each time the equipment is started up. By clicking on the warning, the message will disappear and you could work normally.

7.2 Setting screen

Icons meaning:

7.2.1 Set time and date

Press on the time/date menu. Press successively on the box to select the day, month, year, minute and second adjustment.

Adjust the values by pressing on the icons Increase and Decrease.

7.2.2 Set parameters of temperature control and «Autotuning»

With a value of the appropriate P, I, D parameters, the oven will reach the temperature in the shortest possible time and without exceeding the set temperature. However, those values can be different at high and low temperatures.

The oven comes with the values which are appropriate for most jobs.

These values can be modified, but mathematical knowledge is needed to calculate them. It is much more advisable to perform an "Autotuning" when you need more accuracy.

"Autotuning" automatically adjusts the optimal parameters for a given temperature and load. This process is long because it performs several heating / cooling phases to know the functioning of the environment.

Press on the Control menu. Press successively on the box to select the P, I, D values adjustment and modify the value by pressing on the Increase and Decrease icons.

Attention: It is not recommended to change these parameters. Use the «Autotuning» option to automatically adjust these values.

Select the «Autotuning» temperature with:   (it is recommended to set the most used temperature). There are two calculation methods that are identified with 100 or with 40. The 100 method has a longer duration, but it is more accurate. For most applications, 40 method is more than enough.

Start «Autotuning» by pressing on the icon 

If for some reason you want to stop, press on: 

Note: the «Autotuning» process has a long duration, it is important that the conditions do not vary, therefore do not open the door or handle the oven.


Fig 7.2. Communication

7.2.3 Communication type

The equipment can store readings from time to time and send them to a computer or store them in an USB memory pen drive.

- «PORT» Set the readings storing mode:
 - RS232: Send the temperature readings to a computer by a serial cable (DB9 connector).
 - USB: Store the readings in the internal memory and download later to a USB
- «INTERVAL» Select with ▲ ▼ the register time interval or the data sending.
- «MODE» Select the data to be stored:
 - Start: Save all temperature readings.
 - Set point: Save temperature readings, since the set point temperature has been reached.
- «USB» Update firmware.

Note: The system only saves the last process. To save the process to an USB memory, simply connect the USB to the connector, and data will send automatically to the USB memory. A "sent OK" message will appear. If a process is being performed, this action will be performed after pressing STOP.


Fig 7.2. Temperature alarms

7.2.4 Set the temperature and the temperature alarms

It is recommended to verify calibration of temperature with some regularity (for example, annually) by using for instance a thermometer or a temperature pattern.

If it is necessary to correct the reading as a result of a calibration, one must enter a correction.

To do this, press on the temperature menu.

A box indicating the current correction value appears.

Adjust with keys ▲ ▼ until you get the indicated value on the temperature pattern.

Applied correction will appear in brackets.

7.2.4.1 Temperature alarms

By pressing successively on the temperature correction box, one may enter an over-temperature and a low-temperature limit. Both alarms activate only a message on screen, if the temperature is outside the limits.

7.2.5 Temperature units

Press on the temperature units menu. Press successively to change between Celsius or Farenheit degrees. It is also possible to change the ramp units de °C/min or °C/hour.

7.2.6 Set language

Press on the language menu. Press successively on the box to select among Spanish, French or English.

Adjust by pressing on the Increase and Decrease icons.

7.3 TASK screen

There are two operation modes:

- Working at a fixed temperature.
- Working according to a pre-programmed temperature profile.

Press on  to commute the operation mode (fixed/program)


Fig 7.3.1. Run screen operating at one temperature.

7.3.1 Operation at a fixed temperature

- Select the set point temperature by pressing on the temperature indicator.
- Select a temperature increase ramp if you want to increase it in a slowly way. The value "--" makes that the ramp increases as fast as possible.
- Select an operating time. This time starts to count down when the temperature reaches the set point. The INF value sets a continuous operation.

- Start by pressing on this icon: 

7.3.2 Operation with a temperature profile

- Temperature profiles must be previously programmed from the «PROGRAMMATION» display.
- Select a profile (from 0 to 9) by pressing on keys:  
- Start by pressing on the icon: 
- In the temperature box (top-left corner) the current temperature and the set point are indicated.


Fig 7.3. Run screen operating with temperature profile.

Fig 7.3.2 Working screen data.

Step target temperature	*	Program:	Current step duration:
Current temperature		Nº step on running.	Last ramp run value.
			Total elapsed time of the current step:


- * Point 1 (top) indicates the status of the oven:
 - . if it flashes red, the heating works.
 - . if it is yellow (only in programs) it indicates that the stove is in the resting phase (waiting to reach the heating phase).
- * Point 2 (lower) indicates the data log:
 - . if it flashes green, it indicates that data logging is activated.
 - . if it flashes white, it indicates the frequency of data recording, or what is the same, the recording interval.

7.4 PROGRAMMING screen

Go to PROGRAMMING screen to:

- Set the temperature profiles.
- Set the timed operation for a profile.
(For instance, you can set a time and a week day to start a profile. As start, profile nº 1 on Sunday at 19h00')

7.4.1 Programming a temperature profile:

A temperature profile is a sequence of steps, also called "dells".

To go from the end of a «dell» to the following «dell», it could be done in two ways:

- With a «temperature ramp».
(go from 100°C to 200°C following a ramp at 2°C per minute).
- As fast as possible, "--" indicator (maximum power).

Profile parameters input:

- Select the profile number between 0 and 10 with ▲ ▼ .
- Press on the parameter to be modified. (Temp, ramp, length ...)
- Press on these icons: ▲ ▼ to enter the value.
- Modify another parameter.

Special values:

- If you enter value **OFF** in **temperature** (you have to press Increase key) and the next steps are disabled and not executed.
- If you enter **INF** in **Step length**, the oven remains permanently at this temperature.
- If you enter -- in the **ramp**, the oven makes the temperature rise at the fastest speed possible.

Starting of temperature profile

- To set the profile starting time, press on «—:—» (09:00 in the picture).
- Enter time by pressing ▲ ▼ in **hh:mm** format.
- Press on the days of the week when the profile should be run.
- Accept the programming by pressing on .

Note: If you run a timed program, do not forget to press the RUN key to start the program. During the timeout, it will indicate that it is in Standby mode.

8. USB log

The device has an internal memory for data acquisition (Datalogger).

Once the cycle is finished, it is possible to download this information to an USB memory via an USB connector. The device (USB memory) must be formatted to FAT 32.

Functioning:

- By default, the registrar is not active.

Follow the instructions to change the log interval (see settings>log). Adjust this interval with the conditions of your program.

With the cycle finished:

- Connect the «USB memory» to the USB connector on the front of the device.
- When an «USB memory» is connected, the registration data will be automatically downloaded. A message appears indicating that a connected memory is found and the message "send ok" appears at the end of the download.
- This operation can be performed several times, as long as another cycle has not started (press RUN). When pressing RUN, internal memory of the device empties. The device only stores the last process performed.

Recover data from PC:

- Connect the «USB memory» to a computer to download the data file.
- Important note: The "USB memory" must have a device name. It cannot be blank.

The created file is text type, which can be opened with WordPad, WORD or EXCEL.

9. Maintenance

Before removing the equipment cover, disconnect it from the mains. The control board must only be handled by duly authorized maintenance personnel.

Cleaning:

- Inner enclosure: clean with polish for stainless steel.
- Outer enclosure: clean with a damp cloth or a mild detergent.
- Touch screen: use a cloth and a glass-cleaning liquid.

Annually:

- Contrast the temperature against recognized standart traceable equipments.
- Visually check the status of electrical connections, grounding and connection cable to the mains.

J.P. SELECTA, s.a.u. 2020